

Anti-WNT5A Antibody Picoband®

Catalog Number: A00549-1

About WNT5A

Protein Wnt-5a is a protein that in humans is encoded by the WNT5A gene. This gene is a member of the WNT gene family. The WNTs comprise a large class of secreted proteins that control essential developmental processes such as embryonic patterning, cell growth, migration, and differentiation. It is mapped to 3p14.3. WNT5A is highly expressed in the dermal papilla of depilated skin. This protein plays an essential role in regulating developmental pathways during embryogenesis. And this protein may also play a role in oncogenesis. WNT5A and DVL can promote axon differentiation mediated by the PAR3-PAR6-aPKC complex.

Overview


Product Name	Anti-WNT5A Antibody Picoband®
Reactive Species	Human
Description	Boster Bio Anti-WNT5A Antibody Picoband® catalog # A00549-1. Tested in ELISA, WB applications. This antibody reacts with Human. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P41221

Technical Details

Immunogen	E.coli-derived human WNT5A recombinant protein (Position: D198-D303). Human WNT5A shares 100% amino acid (aa) sequence identity with both mouse and rat WNT5A.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Lyophilized

Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml, Human ELISA, 0.1-0.5 ug/ml, -

Anti-WNT5A Antibody Picoband® (A00549-1) Images



Western blot analysis of WNT5A using anti-WNT5A antibody (A00549-1). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human SiHa whole cell lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-WNT5A antigen affinity purified polyclonal antibody (Catalog # A00549-1) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for WNT5A at approximately 43 kDa. The expected band size for WNT5A is at 42 kDa.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-WNT5A Antibody

For Research Use Only. Not for use in diagnostic procedures.