

## Anti-Factor H CFH Antibody

Catalog Number: A00562-1

### Overview

Product Name	Anti-Factor H CFH Antibody
Reactive Species	Human
Description	Boster Bio Anti-Factor H CFH Antibody catalog # A00562-1. Tested in WB,ICC applications. This antibody reacts with Human.
Application	ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P08603

### Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 ICC: 1:50-1:200

## Anti-Factor H CFH Antibody (A00562-1) Images




Figure 2. Immunocytochemistry staining of CFH using Anti-Factor H CFH Antibody (A00562-1).

ICC staining Factor H in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde




Figure 1. Western blotting validation for Anti-Factor H CFH Antibody A00562-1

Western blot analysis of Factor H on human lung tissue lysate using anti-Factor H antibody at 1/1

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Factor H CFH Antibody