

Anti-Fibronectin FN1 Antibody

Catalog Number: A00564-2


Overview

Product Name	Anti-Fibronectin FN1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Fibronectin FN1 Antibody catalog # A00564-2. Tested in WB,ICC/IF,IHC applications. This antibody reacts with Human.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P02751


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000 ICC: 1:100-1:500 IHC: 1:50-1:200

Anti-Fibronectin FN1 Antibody (A00564-2) Images



ICC staining Fibronectin in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining Fibronectin in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

1 Publications Citing This Product

1. PubMed ID: 33389493, Xu L,Sun N,Li G,Liu L.LncRNA H19 promotes keloid formation through targeting the miR-769-5p/EIF3A pathway.Mol Cell Biochem.2021 Jan 3.doi:10.1007/s11010-020-04024-x.Epub ahead of print.PMID:33389493.

Visit bosterbio.com/anti-fibronectin-fn1-antibody-a00564-2-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fibronectin FN1 Antibody

For Research Use Only. Not for use in diagnostic procedures.