

## Anti-ULK1 (S317) Antibody

Catalog Number: A00584S317

### Overview

Product Name	Anti-ULK1 (S317) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ULK1 (S317) Antibody catalog # A00584S317. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75385

### Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human VPS34.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-ULK1 (S317) Antibody (A00584S317) Images



Western blot (WB) analysis of ULK1 polyclonal antibody at 1:1000 dilution Lane1:L02 whole cell lysate(20ug) Lane2:Panc1 whole cell lysate(20ug) Lane3:HEK293T whole cell lysate(40ug) Lane4:HCC827 whole cell lysate(40ug) Lane5:The Brain tissue lysate of Rat(40ug) Lane6:The Brain tissue lysate of Mouse(40ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ULK1 (S317) Antibody

For Research Use Only. Not for use in diagnostic procedures.