

Anti-ULK1 (S757) Antibody

Catalog Number: A00584S757


Overview

Product Name	Anti-ULK1 (S757) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ULK1 (S757) Antibody catalog # A00584S757. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75385

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 300-350 of Human ULK1.
Isotype	IgG
Form	Liquid
Concentration	1mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-ULK1 (S757) Antibody (A00584S757) Images



Western blot (WB) analysis of ULK1 pAb at 1:500 dilution
Lane1:A549 whole cell lysate(40ug) Lane2:HEK293 whole cell lysate(40ug) Lane3:AML-12 whole cell lysate(40ug)
Lane4:C6 whole cell lysate(30ug) Lane5:Panc1 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ULK1 (S757) Antibody

For Research Use Only. Not for use in diagnostic procedures.