

Anti-ADAMTS13 Antibody

Catalog Number: A00586


Overview

Product Name	Anti-ADAMTS13 Antibody
Reactive Species	Human
Description	Boster Bio Anti-ADAMTS13 Antibody catalog # A00586. Tested in WB,IHC,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q76LX8


Technical Details

Immunogen	Recombinant protein within human ADAMTS13 aa 150-350.
Isotype	IgG
Form	Liquid
Concentration	1mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-ADAMTS13 Antibody (A00586) Images



Western blot analysis of ADAMTS13 on SiHa cell (1), PC-3M cell (2) and A549 cell lysates using anti-ADAMTS13 antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded human liver tissue using anti-ADAMTS13 antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ADAMTS13 Antibody

For Research Use Only. Not for use in diagnostic procedures.