

Anti-Integrin Beta 3 (S778) ITGB3 Antibody

Catalog Number: A00587S778


Overview

Product Name	Anti-Integrin Beta 3 (S778) ITGB3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Integrin Beta 3 (S778) ITGB3 Antibody catalog # A00587S778. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P05106


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 751-800 of Human Integrin beta3.
Isotype	IgG
Form	Liquid
Concentration	1mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-Integrin Beta 3 (S778) ITGB3 Antibody (A00587S778) Images



Western blot (WB) analysis of Integrin beta3 (S778) pAb at 1:500 dilution Lane1:BV2 whole cell lysate(40ug) Lane2:RAW264.7 whole cell lysate(40ug) Lane3:The Uterus tissue lysate of Rat(40ug) Lane4:THP-1 whole cell lysate(40ug) Lane5:A549 whole cell lysate(40ug) Lane6:PC3 whole cell lysate(40ug) Lane7:U-87MG whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of Integrin beta3 (S778) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin Beta 3 (S778) ITGB3 Antibody

For Research Use Only. Not for use in diagnostic procedures.