

## Anti-GATA3 (Phospho-S308) Antibody

Catalog Number: A00593S308


### Overview

Product Name	Anti-GATA3 (Phospho-S308) Antibody
Reactive Species	Human
Description	Boster Bio Anti-GATA3 (Phospho-S308) Antibody catalog # A00593S308. Tested in WB,IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P23771


### Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser308 of human GATA3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 IHC: 1:50-1:200

## Anti-GATA3 (Phospho-S308) Antibody (A00593S308) Images



Western blot analysis of Phospho-GATA3(S308) on human skin lysates using anti-Phospho-GATA3(S308) antibody at 1/1,000 dilution.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using anti-Phospho-GATA3(S308) antibody. Counter stained with hematoxylin.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GATA3 (Phospho-S308) Antibody

For Research Use Only. Not for use in diagnostic procedures.