

## Anti-CD26 DPP4 Antibody

Catalog Number: A00597-1


### Overview

Product Name	Anti-CD26 DPP4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CD26 DPP4 Antibody catalog # A00597-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P27487


### Technical Details

Immunogen	Recombinant fusion protein of human CD26/DPP4(NP_001926.2).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200

## Anti-CD26 DPP4 Antibody (A00597-1) Images



Western blot analysis of extracts of various cell lines, using CD26/DPP4 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.



Immunohistochemistry of paraffin-embedded rat spleen using CD26/DPP4 Rabbit pAb at dilution of 1:25. Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD26 DPP4 Antibody

For Research Use Only. Not for use in diagnostic procedures.