

Anti-Connexin 43 (phospho-S367) GJA1 Antibody

Catalog Number: A00599S367


Overview

Product Name	Anti-Connexin 43 (phospho-S367) GJA1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Connexin 43 (phospho-S367) GJA1 Antibody catalog # A00599S367. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P17302

Technical Details

Immunogen	Synthetic phosphopeptide derived from human Connexin 43 around the phosphorylation site of Serine 367.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-Connexin 43 (phospho-S367) GJA1 Antibody (A00599S367) Images



Western blot (WB) analysis of p-Connexin 43 (S367) pAb at 1:500 dilution Lane1:SGC7901 whole cell lysate(40ug) Lane2:SGC7901 treated with PMA(50ng/ml,10 minutes) whole cell lysate(40ug) Lane3:The Heart tissue lysate of Mouse(40ug) Lane4:The Testis tissue lysate of Rat(40ug)



Immunohistochemistry (IHC) analyzes of p-Connexin 43 (S367) pAb in paraffin-embedded human brain tissue.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Connexin 43 (phospho-S367) GJA1 Antibody

For Research Use Only. Not for use in diagnostic procedures.