

## Anti-HMGCR Antibody

Catalog Number: A00643-2

### Overview


|                      |  |
|----------------------|--|
| Product Name         | Anti-HMGCR Antibody  |
| Reactive Species     | Human, Mouse, Rat  |
| Description          | Boster Bio Anti-HMGCR Antibody catalog # A00643-2. Tested in WB,IP applications. This antibody reacts with Human,Mouse,Rat.                |
| Application          | IP, WB   |
| Clonality            | Polyclonal   |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2   |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit   |
| Uniprot ID           | P04035   |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | recombinant protein                      |
| Cross Reactivity    | No cross reactivity with other proteins. |
| Isotype             | IgG                                      |
| Form                | Liquid                                   |
| Concentration       | 1 mg/ml                                  |
| Purification        | ProA affinity purified                   |
| Suggested Dilutions | WB: 1:1,000-1:2,000                      |

## Anti-HMGCR Antibody (A00643-2) Images

---



Western blot analysis of HMGCR on HeLa cells lysates using anti-HMGCR antibody at 1/1,000 dilution.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-HMGCR Antibody

For Research Use Only. Not for use in diagnostic procedures.