

## Anti-Transcription factor EB TFEB Antibody

Catalog Number: A00662-1


### Overview

Product Name	Anti-Transcription factor EB TFEB Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Transcription factor EB TFEB Antibody catalog # A00662-1. Tested in ELISA, IHC, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P19484


### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human TFEB. AA range:10-59
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:20000 IF 1:50-200

## Anti-Transcription factor EB TFEB Antibody (A00662-1) Images



Western blotting validation for Anti-Transcription factor EB TFEB Antibody A00662-1 Western Blot (WB) analysis of Mouse Lung Mouse Muscle cells using TFEB polyclonal antibody. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



Western blotting validation for Anti-Transcription factor EB TFEB Antibody A00662-1 Western Blot (WB) analysis of Mouse Kidney lysis using TFEB antibody. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Transcription factor EB TFEB Antibody

For Research Use Only. Not for use in diagnostic procedures.