

Anti-CXCR5 Antibody

Catalog Number: A00663

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-CXCR5 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CXCR5 Antibody catalog # A00663. Tested in WB,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P32302

Technical Details


Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200 FC: 1:50-1:100

--	--

Anti-CXCR5 Antibody (A00663) Images



Western blot analysis of CXCR5 on Daudi cell using anti-CXCR5 antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded rat kidney tissue using anti-CXCR5 antibody. Counter stained with hematoxylin.

[Anti-CXCR5 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.