

Anti-HPRT HPRT1 Antibody

Catalog Number: A00668


Overview

Product Name	Anti-HPRT HPRT1 Antibody
Reactive Species	Human, Mouse, Rat, Zebrafish
Description	Boster Bio Anti-HPRT HPRT1 Antibody catalog # A00668. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat,Zebrafish.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P00492


Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200 IP: 1:10-1:50

Anti-HPRT HPRT1 Antibody (A00668) Images



Western blot analysis of HPRT on different lysates using anti-HPRT antibody at 1/500 dilution. Positive control: Lane 1: HeLa Lane 2: A431 Lane 3: MCF-7 Lane 4: 293 Lane 5: Rat kidney Lane 6: Rat brain Lane 7: Mouse testis Lane 8: Mouse colon



Immunohistochemical analysis of paraffin-embedded rat brain tissue using anti-HPRT antibody. Counter stained with hematoxylin.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HPRT HPRT1 Antibody

For Research Use Only. Not for use in diagnostic procedures.