

Anti-p68 RNA Helicase (phospho-Y593) DDX5 Antibody

Catalog Number: A00670Y593


Overview

Product Name	Anti-p68 RNA Helicase (phospho-Y593) DDX5 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-p68 RNA Helicase (phospho-Y593) DDX5 Antibody catalog # A00670Y593. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P17844

Technical Details

Immunogen	Synthetic peptide derived from human p68 RNA Helicase around the phosphorylation site of Y593.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-p68 RNA Helicase (phospho-Y593) DDX5 Antibody (A00670Y593) Images



Western blot (WB) analysis of p68 RNA Helicase (phospho-Y593) polyclonal antibody at 1:500 dilution Lane1:MCF-7 whole cell lysate(40µg) Lane2:A375 whole cell lysate(40µg) Lane3:U-87MG whole cell lysate(40µg) Lane4:CT-26 whole cell lysate(40µg) Lane5:C6 whole cell lysate(40µg)

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p68 RNA Helicase (phospho-Y593) DDX5 Antibody

For Research Use Only. Not for use in diagnostic procedures.