

Anti-HSP27 (phospho-S78) HSPB1 Antibody

Catalog Number: A00676S78


Overview

Product Name	Anti-HSP27 (phospho-S78) HSPB1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HSP27 (phospho-S78) HSPB1 Antibody catalog # A00676S78. Tested in IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P04792

Technical Details

Immunogen	Synthetic phosphopeptide derived from human HSP27 around the phosphorylation site of Serine 78.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-HSP27 (phospho-S78) HSPB1 Antibody (A00676S78) Images



Immunohistochemistry (IHC) analysis of p-HSP27 (Ser78) pAb in paraffin-embedded human breast carcinoma tissue

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HSP27 (phospho-S78) HSPB1 Antibody

For Research Use Only. Not for use in diagnostic procedures.