

## Anti-Glucagon GCG Antibody

Catalog Number: A00678-1


### Overview

Product Name	Anti-Glucagon GCG Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Glucagon GCG Antibody catalog # A00678-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01275


### Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200

## Anti-Glucagon GCG Antibody (A00678-1) Images



Immunohistochemical analysis of paraffin-embedded human pancreas tissue using anti-Glucagon antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded mouse pancreas tissue using anti-Glucagon antibody. Counter stained with hematoxylin.

## 1 Publications Citing This Product

1. PubMed ID: 32328034, Huo K,Li X,Hu W,Song X,Zhang D,Zhang X,Chen X,Yuan J,Zuo J,Wang X.RFRP-3, the Mammalian Ortholog of GnIH, Is a Novel Modulator Involved in Food Intake and Glucose Homeostasis.Front Endocrinol (Lausanne). 2020 Apr 9;11:194.doi:10.3389/fendo.2020.00194.PMID

Visit [bosterbio.com/anti-glucagon-gcg-antibody-a00678-1-boster.html](https://bosterbio.com/anti-glucagon-gcg-antibody-a00678-1-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Glucagon GCG Antibody

For Research Use Only. Not for use in diagnostic procedures.