

Anti-Mucin-5B MUC5B Antibody

Catalog Number: A00719

Overview

Product Name	Anti-Mucin-5B MUC5B Antibody
Reactive Species	Human
Description	Boster Bio Anti-Mucin-5B MUC5B Antibody catalog # A00719. Tested in IHC, IF applications. This antibody reacts with Human.
Application	IF, IHC
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9HC84

Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 1810-1890
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	MUC5B antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:50-300 IF 1:50-200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Mucin-5B MUC5B Antibody

For Research Use Only. Not for use in diagnostic procedures.