

Anti-Nibrin (phospho-S343) NBN Antibody

Catalog Number: A00732S343


Overview

Product Name	Anti-Nibrin (phospho-S343) NBN Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Nibrin (phospho-S343) NBN Antibody catalog # A00732S343. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O60934


Technical Details

Immunogen	Synthetic phosphopeptide derived from human Nibrin around the phosphorylation site of Serine 343.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-Nibrin (phospho-S343) NBN Antibody (A00732S343) Images



Western blot (WB) analysis of p-Nibrin (S343) pAb at 1:500 dilution Lane1:HEK293T whole cell lysate(40ug) Lane2:HEK293T treated with UV for 5 minutes then repair for 1 hours whole cell lysate(40ug) Lane3:HEK293T treated with UV for 5 minutes then repair for 3 hours whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of p-p95/NBS1 (S343) pAb in paraffin-embedded human tonsil tissue.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Nibrin (phospho-S343) NBN Antibody

For Research Use Only. Not for use in diagnostic procedures.