

## Anti-p95/NBS1 (Phospho-S343) NBN Antibody

Catalog Number: A00732S343-1


### Overview

Product Name	Anti-p95/NBS1 (Phospho-S343) NBN Antibody
Reactive Species	Human
Description	Boster Bio Anti-p95/NBS1 (Phospho-S343) NBN Antibody catalog # A00732S343-1. Tested in WB,ICC/IF,IHC,IP applications. This antibody reacts with Human.
Application	IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O60934


### Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser343 of human p95/NBS1.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 ICC: 1:50-1:200 IHC: 1:50-1:200

## Anti-p95/NBS1 (Phospho-S343) NBN Antibody (A00732S343-1) Images



ICC staining Phospho-p95/NBS1(S343) in PC-3M cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



Immunohistochemical analysis of paraffin-embedded mouse testis tissue using anti-Phospho-p95/NBS1(S343) antibody. Counter stained with hematoxylin.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p95/NBS1 (Phospho-S343) NBN Antibody

For Research Use Only. Not for use in diagnostic procedures.