

## Anti-Stat4 (phospho-Y693) Antibody

Catalog Number: A00734Y693


### Overview

Product Name	Anti-Stat4 (phospho-Y693) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Stat4 (phospho-Y693) Antibody catalog # A00734Y693. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14765


### Technical Details

Immunogen	Synthetic phosphopeptide derived from human Stat4 around the phosphorylation site of Tyrosine 693.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IP: 1:50-1:200

## Anti-Stat4 (phospho-Y693) Antibody (A00734Y693) Images



Western blot (WB) analysis of p-Stat4 (Y693) pAb at 1:500 dilution Lane1:MCF-7 whole cell lysate(40ug) Lane2:MCF-7 treated with IFN-gamma(100ng/ml,15 minutes) whole cell lysate(40ug) Lane3:MCF-7 treated with IFN-gamma(100ng/ml,30 minutes) whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of p-STAT4 (Y693) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Stat4 (phospho-Y693) Antibody

For Research Use Only. Not for use in diagnostic procedures.