

## Anti-TNF-beta LTA Antibody

Catalog Number: A00746-2

### Overview

Product Name	Anti-TNF-beta LTA Antibody
Reactive Species	Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-TNF-beta LTA Antibody catalog # A00746-2. Tested in ELISA, WB applications. This antibody reacts with Human, Monkey, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01374


### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human TNF-beta. AA range:17-66
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:40000

## Anti-TNF-beta LTA Antibody (A00746-2) Images



Western Blot analysis of various cells using TNF-beta Polyclonal Antibody diluted at 1:1000. Secondary antibody was diluted at 1:20000



Western blot analysis of lysate from COS-7cell, using TNF-beta antibodies.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNF-beta LTA Antibody

For Research Use Only. Not for use in diagnostic procedures.