

Anti-ZAP70 (T286) Antibody

Catalog Number: A00754T286


Overview

Product Name	Anti-ZAP70 (T286) Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-ZAP70 (T286) Antibody catalog # A00754T286. Tested in WB,IHC applications. This antibody reacts with Human,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P43403

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 250-300 of Human ZAP70.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-ZAP70 (T286) Antibody (A00754T286) Images



Western blot (WB) analysis of ZAP70 (T286) at 1:500 dilution
Lane1:THP-1 whole cell lysate(40ug) Lane2:Myla2059 whole cell lysate(40ug) Lane3:HuT78 whole cell lysate(40ug)
Lane4:The Thymus tissue lysate of Rat(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ZAP70 (T286) Antibody

For Research Use Only. Not for use in diagnostic procedures.