

## Anti-ZAP70 (Tyr292) Antibody

Catalog Number: A00754Y292


### Overview

Product Name	Anti-ZAP70 (Tyr292) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-ZAP70 (Tyr292) Antibody catalog # A00754Y292. Tested in WB applications. This antibody reacts with Human,Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P43403

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 280-320 of Human ZAP70.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-ZAP70 (Tyr292) Antibody (A00754Y292) Images



Lane1: Hela whole cell lysate  
Lane2:Hep3B whole cell lysate  
Lane3: sp2/O whole cell lysate  
Lane4:Raw264.7 whole cell lysate  
ZAP-70(Tyr292) pAb at 1:500 dilution

Western blot (WB) analysis of ZAP70 (Tyr292) pAb

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-ZAP70 (Tyr292) Antibody

For Research Use Only. Not for use in diagnostic procedures.