

Anti-Phospho-ZAP-70 (Y493) Antibody

Catalog Number: A00754Y493

About ZAP70

Couples tyrosine kinase signals with the activation of the Rho/Rac GTPases, thus leading to cell differentiation and/or proliferation.

Ana V. Miletic, J. Biol. Chem., Dec 2006; 281: 38257 - 38265.

Christopher C. Stebbins, Mol. Cell. Biol., Sep 2003; 23: 6291 - 6299.

Vered Schapira, Cancer Res., Jun 2006; 66: 6183 - 6191.

Maya Groysman, J. Biol. Chem., Dec 2002; 277: 50121 - 50130.

Overview

Product Name	Anti-Phospho-ZAP-70 (Y493) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-ZAP-70 (Y493) Antibody catalog # A00754Y493. Tested in ELISA, IHC, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P43403

Technical Details

Immunogen	Synthesized peptide derived from human ZAP-70 around the phosphorylation site of Y493.
Predicted Reactive Species	Boar, Bovine, Canine, Golden Hamster
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-

	specific immunogen.
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:5000</p>

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-ZAP-70 (Y493) Antibody