

Anti-CCR4 Antibody Picoband® FITC Conjugated

Catalog Number: A00755-FITC

About CCR4

CCR4, C-C chemokine receptor type 4 is a protein that in humans is encoded by the CCR4 gene. CCR4 has also recently been designated CD194 (cluster of differentiation 194). The protein encoded by this gene belongs to the G protein-coupled receptor family. It is a receptor for the following CC chemokines. The CCR4 gene is mapped to chromosome 3. Chemokines are a group of small structurally related proteins that regulate cell trafficking of various types of leukocytes. The chemokines also play fundamental roles in the development, homeostasis, and function of the immune system, and they have effects on cells of the central nervous system as well as on endothelial cells involved in angiogenesis or angiostasis.

Overview

| | |
|----------------------|--|
| Product Name | Anti-CCR4 Antibody Picoband® FITC Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P51679 |

Technical Details

| | |
|---------------------|---|
| Immunogen | A synthetic peptide corresponding to a sequence at the C-terminus of human CCR4, which shares 93.3% amino acid (aa) sequence identity with both mouse and rat CCR4. |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users. |

1 Publications Citing This Product

1. PubMed ID: 24416340, Zhao C, Bu X, Wang W, Ma T, Ma H. Plos One. 2014 Jan 8;9(1):E85058. Doi: 10.1371/Journal.Pone.0085058. Ecollection 2014. Gec-Derived Sfrp5 Inhibits Wnt5A-Induced Macrophage Chemotaxis And Activation.

Visit bosterbio.com/anti-ccr4-picoband-trade-antibody-a00755-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CCR4 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.