

## Anti-TGF Beta RII (I122) TGFBR2 Antibody

Catalog Number: A00759


### Overview

Product Name	Anti-TGF Beta RII (I122) TGFBR2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TGF Beta RII (I122) TGFBR2 Antibody catalog # A00759. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P37173

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 100-150 of Human TGFbeta RII.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-TGF Beta RII (I122) TGFBR2 Antibody (A00759) Images



Western blot (WB) analysis of TGFbeta RII (I122) polyclonal antibody at 1:500 dilution Lane1:Hela whole cell lysate Lane2:Raw264.7 whole cell lysate Lane3:H9C2 whole cell lysate

### 1 Publications Citing This Product

1. PubMed ID: -, xiangrong chen, Yile Zeng, Fan Wang et al. ADAM17 Aggravates the Inflammatory Response by Modulating Microglia Polarization Through the TGF-beta1/Smad Pathway Following Experimental Traumatic Brain Injury, 04 February 2021, PREPRINT (Version 1) available at

Visit [bosterbio.com/anti-tgf-beta-rii-i122-tgfr2-antibody-a00759-boster.html](https://bosterbio.com/anti-tgf-beta-rii-i122-tgfr2-antibody-a00759-boster.html) to see all 1 publications.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TGF Beta RII (I122) TGFBR2 Antibody

For Research Use Only. Not for use in diagnostic procedures.