

Anti-GATA-6 (R493) Antibody

Catalog Number: A00778


Overview

Product Name	Anti-GATA-6 (R493) Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-GATA-6 (R493) Antibody catalog # A00778. Tested in WB applications. This antibody reacts with Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92908

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 460-510 of Human GATA-6.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-GATA-6 (R493) Antibody (A00778) Images



Western blot (WB) analysis of GATA-6 (R493) pAb at 1:500 dilution
Lane1: The Heart tissue lysate of Mouse(40ug)
Lane2: The Heart tissue lysate of Rat(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GATA-6 (R493) Antibody

For Research Use Only. Not for use in diagnostic procedures.