

Anti-Fibrinogen Gamma FGG Antibody

Catalog Number: A00790


Overview

Product Name	Anti-Fibrinogen Gamma FGG Antibody
Reactive Species	Human
Description	Boster Bio Anti-Fibrinogen Gamma FGG Antibody catalog # A00790. Tested in WB,IHC,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P02679


Technical Details

Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:5,000-1:10,000 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-Fibrinogen Gamma FGG Antibody (A00790) Images



Western blot analysis of Fibrinogen gamma chain on human plasma using anti-Fibrinogen gamma chain at 1/5,000 dilution.



Immunohistochemical analysis of paraffin-embedded human colon tissue using anti-Fibrinogen gamma chain antibody. Counter stained with hematoxylin.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fibrinogen Gamma FGG Antibody

For Research Use Only. Not for use in diagnostic procedures.