

Anti-GSK3 Beta (Phospho-Ser9) GSK3B Antibody

Catalog Number: A00791S9-1


Overview

Product Name	Anti-GSK3 Beta (Phospho-Ser9) GSK3B Antibody
Reactive Species	Human
Description	Boster Bio Anti-GSK3 Beta (Phospho-Ser9) GSK3B Antibody catalog # A00791S9-1. Tested in WB,ICC/IF,IHC applications. This antibody reacts with Human.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P49841


Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser9 of human GSK3 beta.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 ICC: 1:50-1:200 IHC: 1:50-1:200

Anti-GSK3 Beta (Phospho-Ser9) GSK3B Antibody (A00791S9-1) Images



Western blot analysis of Phospho-GSK3 beta(Ser 9) on different lysates using anti-Phospho-GSK3 beta(Ser 9) antibody at 1/1,000 dilution. Positive control: Lane 1: HeLa
Lane 2: MCF-7



ICC staining Phospho-GSK3 beta(Ser 9) in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GSK3 Beta (Phospho-Ser9) GSK3B Antibody

For Research Use Only. Not for use in diagnostic procedures.