

Anti-DMPK (R30) Antibody

Catalog Number: A00797-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-DMPK (R30) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DMPK (R30) Antibody catalog # A00797-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q09013

Technical Details

Immunogen	Synthetic peptide, corresponding to the N-terminus of Human DMPK.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

--	--

Anti-DMPK (R30) Antibody (A00797-1) Images



Western blot (WB) analysis of DMPK (R30) pAb at 1:500 dilution Lane1:CT26 whole cell lysate(40ug) Lane2:The Brain tissue lysate of Rat(40ug) Lane3:A549 whole cell lysate(40ug) Lane4:L02 whole cell lysate(40ug) Lane5:HEK293T whole cell lysate(20ug)

Anti-DMPK (R30) Antibody

For Research Use Only. Not for use in diagnostic procedures.