

## Anti-DMPK Antibody Picoband®

Catalog Number: A00797-2

### About DMPK

The protein encoded by this gene is a serine/threonine protein kinase that contains coiled-coil and C-terminal membrane association domains. In the embryonic mouse, it is found in cardiac and skeletal myocytes where it appears to play a role in myogenesis. In adults, the transcript is localized to several tissues including brain, heart, and skeletal and smooth muscle, and a function in cytoskeletal remodeling has been described. Transcripts with expanded CUG repeats in the 3' untranslated region mediate alternative splicing of several genes and sequester RNA binding proteins and RNA transcripts that contain CAG repeats, resulting in myotonic dystrophy, an autosomal dominant neuromuscular disorder. Alternative splicing results in multiple protein coding and non-coding transcript variants.

### Overview

Product Name	Anti-DMPK Antibody Picoband®
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-DMPK Antibody Picoband® catalog # A00797-2. Tested in WB applications. This antibody reacts with Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P54265


### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human DMPK.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Purification	Immunogen affinity purified.

Suggested Dilutions

Western blot, 0.25-0.5 ug/ml, Mouse, Rat

## Anti-DMPK Antibody Picoband® (A00797-2) Images



Western blot analysis of DMPK using anti-DMPK antibody (A00797-2). Electrophoresis was performed on a 10% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: rat heart tissue lysates, Lane 2: rat skeletal muscle tissue lysates, Lane 3: mouse heart tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-DMPK antigen affinity purified polyclonal antibody (A00797-2) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1% Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for DMPK at approximately 70 kDa. The expected band size for DMPK is at 69 kDa.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-DMPK Antibody

For Research Use Only. Not for use in diagnostic procedures.