

Anti-TBX5 Monoclonal Antibody

Catalog Number: A00803


Overview

Product Name	Anti-TBX5 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-TBX5 Monoclonal Antibody catalog # A00803. Tested in WB, ELISA applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 7B11
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q99593

Technical Details

Immunogen	Purified recombinant fragment of TBX5 expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Ascitic fluid
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-TBX5 Monoclonal Antibody (A00803) Images



Western Blot analysis using TBX5 Monoclonal Antibody against HepG2 cell lysate (1).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TBX5 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.