

Anti-ATG12 (S129) Antibody

Catalog Number: A00820S129


Overview

Product Name	Anti-ATG12 (S129) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ATG12 (S129) Antibody catalog # A00820S129. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O94817


Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human ATG12.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-ATG12 (S129) Antibody (A00820S129) Images



Western blot (WB) analysis of ATG12 (S129) polyclonal antibody at 1:500 dilution Lane1:C6 whole cell lysate(40ug) Lane2:CT26 whole cell lysate(40ug) Lane3:HepG2 whole cell lysate(40ug) Lane4:HCT116 whole cell lysate(40ug)



Immunohistochemistry of paraffin-embedded Human Liver cancer using ATG12 antibody at dilution of 1:50.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ATG12 (S129) Antibody

For Research Use Only. Not for use in diagnostic procedures.