

Anti-Steroid 21-hydroxylase CYP21A2 Antibody

Catalog Number: A00821-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-Steroid 21-hydroxylase CYP21A2 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Steroid 21-hydroxylase CYP21A2 Antibody catalog # A00821-1. Tested in ELISA, IF, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P08686

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 21A2. AA range:151-200
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml

Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 IF 1:200-1:1000 ELISA 1:20000

Anti-Steroid 21-hydroxylase CYP21A2 Antibody (A00821-1) Images



Western blot analysis of 293 cells using CYP21A2 Polyclonal Antibody diluted at 1:1000

[Anti-Steroid 21-hydroxylase CYP21A2 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.