

Anti-Musashi 1 LRP1 Antibody

Catalog Number: A00829


Overview

| | |
|----------------------|--|
| Product Name | Anti-Musashi 1 LRP1 Antibody |
| Reactive Species | Human |
| Description | Boster Bio Anti-Musashi 1 LRP1 Antibody catalog # A00829. Tested in WB,ICC/IF,IHC,IP,Flow Cytometry applications. This antibody reacts with Human. |
| Application | Flow Cytometry, IP, IF, IHC, ICC, WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q07954 |


Technical Details

| | |
|---------------------|--|
| Immunogen | recombinant protein |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | ProA affinity purified |
| Suggested Dilutions | WB: 1:1000 ICC: 1:50-1:200 IHC: 1:50-1:200 FC: 1:50-1:200 |

Anti-Musashi 1 LRP1 Antibody (A00829) Images



ICC staining Musashi 1 in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining Musashi 1 in A549 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Musashi 1 LRP1 Antibody

For Research Use Only. Not for use in diagnostic procedures.