

Anti-ATF2 (Phospho-T71) LRP1 Antibody

Catalog Number: A00829T71


Overview

Product Name	Anti-ATF2 (Phospho-T71) LRP1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-ATF2 (Phospho-T71) LRP1 Antibody catalog # A00829T71. Tested in WB,ICC/IF,IP applications. This antibody reacts with Human.
Application	IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q07954


Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Thr71 of human ATF2
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000 ICC: 1:50-1:200

Anti-ATF2 (Phospho-T71) LRP1 Antibody (A00829T71) Images



ICC staining phospho-ATF2(T71) in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining phospho-ATF2(T71) in NIH/3T3 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ATF2 (Phospho-T71) LRP1 Antibody

For Research Use Only. Not for use in diagnostic procedures.