

## Anti-Rel B (S552) Antibody

Catalog Number: A00836S552


### Overview

Product Name	Anti-Rel B (S552) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Rel B (S552) Antibody catalog # A00836S552. Tested in WB applications. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q01201

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 530-580 of Human Rel B.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

## Anti-Rel B (S552) Antibody (A00836S552) Images



Western blot analysis of Rel B expression in HepG2 cell lysates. The lane on the left is treated with the antigen-specific peptide.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rel B (S552) Antibody

For Research Use Only. Not for use in diagnostic procedures.