

Anti-Annexin A2 ANXA2 Antibody

Catalog Number: A00868


Overview

Product Name	Anti-Annexin A2 ANXA2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Annexin A2 ANXA2 Antibody catalog # A00868. Tested in WB,IHC,ICC/IF,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P07355


Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200 ICC: 1:100-1:200 FC: 1:50-1:200

Anti-Annexin A2 ANXA2 Antibody (A00868) Images



Western blot analysis of Annexin A2 on different cells lysates using anti-Annexin A2 antibody at 1/500 dilution. Positive control. Lane 1: Hela Lane 2: Human kidney Lane 3: Murine testis



ICC staining Annexin A2 in Hela cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Annexin A2 ANXA2 Antibody

For Research Use Only. Not for use in diagnostic procedures.