

Anti-IKK Gamma (H81) IKBKG Antibody

Catalog Number: A00874-3


Overview

Product Name	Anti-IKK Gamma (H81) IKBKG Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IKK Gamma (H81) IKBKG Antibody catalog # A00874-3. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y6K9


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 50-100 of Human IKKgamma.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

Anti-IKK Gamma (H81) IKBKG Antibody (A00874-3) Images



Western blot (WB) analysis of IKBKG pAb at 1:500 dilution
Lane1:MEF whole cell lysate(40ug) Lane2:H9C2 whole cell lysate(40ug) Lane3:H1792 whole cell lysate(40ug)
Lane4:HEK293T whole cell lysate(40ug) Lane5:Panc1 whole cell lysate (40ug)



Immunohistochemistry (IHC) analyzes of IKKgamma (H81) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IKK Gamma (H81) IKBKG Antibody

For Research Use Only. Not for use in diagnostic procedures.