

Anti-Gfi-1 Monoclonal Antibody

Catalog Number: A00888


Overview

Product Name	Anti-Gfi-1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Gfi-1 Monoclonal Antibody catalog # A00888. Tested in WB, ELISA applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 5D7
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q99684

Technical Details

Immunogen	Purified recombinant fragment of human Gfi-1 expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Ascitic fluid
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-Gfi-1 Monoclonal Antibody (A00888) Images



Western Blot analysis using Gfi-1 Monoclonal Antibody against HEK293 (1) and GFI1-hlgGfC transfected HEK293 (2) cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Gfi-1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.