

Anti-Steroidogenic factor 1 NR5A1 Antibody

Catalog Number: A00891


Overview

Product Name	Anti-Steroidogenic factor 1 NR5A1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Steroidogenic factor 1 NR5A1 Antibody catalog # A00891. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13285


Technical Details

Immunogen	Recombinant fusion protein of human Steroidogenic factor-1 (SF-1/Steroidogenic factor-1 (SF-1/NR5A1))(NP_004950.2).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

Anti-Steroidogenic factor 1 NR5A1 Antibody (A00891) Images



Western blot analysis of extracts of Mouse testis, using Steroidogenic factor-1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.



Western blot analysis of extracts of various cell lines, using Steroidogenic factor-1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Steroidogenic factor 1 NR5A1 Antibody

For Research Use Only. Not for use in diagnostic procedures.