

Anti-Protein lin-28 homolog B LIN28B Antibody

Catalog Number: A00896


Overview

Product Name	Anti-Protein lin-28 homolog B LIN28B Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Protein lin-28 homolog B LIN28B Antibody catalog # A00896. Tested in WB,ICC/IF applications. This antibody reacts with Human,Mouse.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6ZN17

Technical Details

Immunogen	Recombinant fusion protein of human LIN28B(NP_001004317.1).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ICC/IF: 1:50-1:200

Anti-Protein lin-28 homolog B LIN28B Antibody (A00896) Images



Western blot analysis of extracts of mouse liver, using LIN28B antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Protein lin-28 homolog B LIN28B Antibody

For Research Use Only. Not for use in diagnostic procedures.