

Anti-Epoxyde hydrolase 1 EPHX1 Antibody

Catalog Number: A00899


Overview

Product Name	Anti-Epoxyde hydrolase 1 EPHX1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Epoxyde hydrolase 1 EPHX1 Antibody catalog # A00899. Tested in WB applications. This antibody reacts with Human,Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P07099

Technical Details

Immunogen	Recombinant fusion protein of human EPHX1(NP_000111.1).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

Anti-Epoxyde hydrolase 1 EPHX1 Antibody (A00899) Images



Western blot analysis of extracts of various cell lines, using EPHX1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Epoxyde hydrolase 1 EPHX1 Antibody

For Research Use Only. Not for use in diagnostic procedures.