

## Anti-Tenascin-C TNC Antibody

Catalog Number: A00936


### Overview

Product Name	Anti-Tenascin-C TNC Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Tenascin-C TNC Antibody catalog # A00936. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P24821

### Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human Tenascin-C.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-Tenascin-C TNC Antibody (A00936) Images



Western blot (WB) analysis of Tenascin-C polyclonal antibody at 1:500 dilution Lane1:MCF-7 cell lysate Lane2:NIH-3T3 cell lysate Lane3:PC12 cell lysate

## 1 Publications Citing This Product

1. PubMed ID: 32648862, Cai J,Xie X,Li D,Wang L,Jiang J,Mo X,Zhao J.A novel knitted scaffold made of microfiber/nanofiber core-sheath yarns for tendon tissue engineering.Biomater Sci.2020 Aug 11;8(16):4413-4425.doi:10.1039/ d0bm00816h.PMID:32648862.

Visit [bosterbio.com/anti-tenascin-c-tnc-antibody-a00936-boster.html](https://bosterbio.com/anti-tenascin-c-tnc-antibody-a00936-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Tenascin-C TNC Antibody

For Research Use Only. Not for use in diagnostic procedures.