

Anti-Interferon regulatory factor 5 IRF5 Antibody

Catalog Number: A00958


Overview

| | |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Anti-Interferon regulatory factor 5 IRF5 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-Interferon regulatory factor 5 IRF5 Antibody catalog # A00958. Tested in WB,ICC/IF,IHC,IP,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat. |
| Application | Flow Cytometry, IP, IF, IHC, ICC, WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q13568 |


Technical Details

| | |
|---------------------|-----------------------------------------------------------------------------|
| Immunogen | recombinant protein |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | ProA affinity purified |
| Suggested Dilutions | WB: 1:1,000-1:5,000 ICC: 1:50-1:200 IHC: 1:50-1:200 FC: 1:50-1:100 |

Anti-Interferon regulatory factor 5 IRF5 Antibody (A00958) Images



Western blot analysis of IRF5 on THP-1 cells lysates using anti-IRF5 antibody at 1/1,000 dilution.



ICC staining IRF5 in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Interferon regulatory factor 5 IRF5 Antibody

For Research Use Only. Not for use in diagnostic procedures.