

## Anti-p130 Cas (P160) BCAR1 Antibody

Catalog Number: A00960-1

### Overview


Product Name	Anti-p130 Cas (P160) BCAR1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-p130 Cas (P160) BCAR1 Antibody catalog # A00960-1. Tested in IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P56945

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 140-180 of Human p130 Cas.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

## Anti-p130 Cas (P160) BCAR1 Antibody (A00960-1) Images

---



Immunohistochemistry (IHC) analyzes of p130 Cas (P160) pAb in paraffin-embedded human breast carcinoma tissue.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

---

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p130 Cas (P160) BCAR1 Antibody

For Research Use Only. Not for use in diagnostic procedures.