

Anti-p130 Cas (phospho-Y249) BCAR1 Antibody

Catalog Number: A00960Y249-1


Overview

Product Name	Anti-p130 Cas (phospho-Y249) BCAR1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-p130 Cas (phospho-Y249) BCAR1 Antibody catalog # A00960Y249-1. Tested in WB,IHC,ICC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P56945


Technical Details

Immunogen	Synthetic phosphopeptide derived from human p130 Cas around the phosphorylation site of Tyrosine 249.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 ICC: 1:50-1:200

Anti-p130 Cas (phospho-Y249) BCAR1 Antibody (A00960Y249-1) Images



Western blot (WB) analysis of p-p130 Cas (Y249) polyclonal antibody at 1:500 dilution Lane1:HEK293T cell lysate treated with EGF(0.1ng/ML,30mins) Lane2:Raw264.7 cell lysate treated with EGF(0.1ng/ML,30mins) Lane3:PC12 cell lysate treated with EGF(0.1ng/ML,30mins)



Immunohistochemistry (IHC) analyzes of p-p130 Cas (Y249) pAb in paraffin-embedded human breast carcinoma tissue at 1:50, showing cytoplasmic staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p130 Cas (phospho-Y249) BCAR1 Antibody

For Research Use Only. Not for use in diagnostic procedures.