

Anti-D4DR DRD4 Antibody

Catalog Number: A00998-1


Overview

Product Name	Anti-D4DR DRD4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-D4DR DRD4 Antibody catalog # A00998-1. Tested in WB, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P21917


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human DRD4. AA range:355-404
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:40000


Anti-D4DR DRD4 Antibody (A00998-1) Images




Immunofluorescence analysis of MCF7 cells, using DRD4 Antibody. The picture on the right is blocked with the synthesized peptide.




Western Blot analysis of various cells using D4DR Polyclonal Antibody



Western Blot analysis of K562 cells using D4DR Polyclonal Antibody



Western blot analysis of lysates from MCF-7 and HUVEC cells, using DRD4 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from RAW264.7 cells using DRD4 antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-D4DR DRD4 Antibody

For Research Use Only. Not for use in diagnostic procedures.