

## Anti-RAGE (Q286) MOK Antibody

Catalog Number: A01004


### Overview

Product Name	Anti-RAGE (Q286) MOK Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-RAGE (Q286) MOK Antibody catalog # A01004. Tested in IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UQ07

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 254-308 of Human RAGE.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

## Anti-RAGE (Q286) MOK Antibody (A01004) Images



Immunohistochemistry (IHC) analyzes of RAGE (Q286) pAb in paraffin-embedded human tonsil cancer tissue at 1:50.

## 3 Publications Citing This Product

1. PubMed ID: 10.1007/s12011-021-02969-8, Lead Exposure in Developmental Ages Promotes Abeta Accumulation by Disturbing Abeta Transportation in Blood-Cerebrospinal Fluid Barrier/Blood-Brain Barriers and Impairing Abeta Clearance in the Liver
2. PubMed ID: 10.1016/j.toxlet.2008.03.014, Proteomic alteration in lung tissue of rats exposed to cigarette smoke
3. PubMed ID: 10.1016/j.jcma.2018.03.013, Involvement of S100A4/Mts1 and associated proteins in the protective effect of fluoxetine against MCT - Induced pulmonary hypertension in rats

Visit [bosterbio.com/anti-rage-q286-mok-antibody-a01004-boster.html](https://bosterbio.com/anti-rage-q286-mok-antibody-a01004-boster.html) to see all 3 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RAGE (Q286) MOK Antibody

For Research Use Only. Not for use in diagnostic procedures.